



Siron® LS

Coagulation test for determination of activated Partial Thromboplastin Time (aPTT). The fluid reagent consists of a watery suspension of phospholipids and ellagic acid. Siron® is sensitive to heparin.

Test principle

The sample is incubated with Siron® LS. At 37°C the reaction with 0,025 molar CaCl_2 is started and the time until formation of a clot measured.

Reagents and stability (after opening)

Siron® LS

14 days RT
3 weeks at 12°C
2 months at 2-8°C

Analyzer applications:	on request
-------------------------------	-------------------

Package size:	2 x 4 ml
----------------------	-----------------

Article no.:	5035105
---------------------	----------------

Package size:	10 x 4 ml
----------------------	------------------

Article no.:	5035107
---------------------	----------------

Package size:	10 x 10 ml
----------------------	-------------------

Article no.:	5035109
---------------------	----------------